

Underivatized Sympathomimetic Amines (Basic Drugs) Rtx[®]-5 Amine

30m, 0.25mm ID, 0.50 μ m Rtx[®]-5 Amine (cat.# 12338)

split mode, split vent flow rate 45mL/min.

Oven temp.: 100°C to 310°C @ 10°C/min.

Varian 3400 GC coupled with Varian Saturn 2000 MS detector

Mass spectroscopy data collected using a scan range of 40 amu through 600 amu.

Ionization performed in the EI Auto mode.

